INTENDED USE

SepaSperm® Solution is a ready to use medium used for separation of motile sperm from seminal fluid with density gradient method.

COMPOSITION:

SepaSperm® Solution medium contains:

- Silane-coated silica particles
- HEPES-buffered
- EBSS (Earle's balanced salt solution)

GENERAL INFORMATION:

SepaSperm® Solution are available in bottles of 20ml, 50 ml and 100 ml.

INSTRUCTIONS FOR USE:

CAUTION: The detailed usage of SepaSperm® Solution is determined based on each laboratory's procedure and protocols that are optimized for individual medical purposes.

Preparation

• Allow the semen to liquefy at 37°C and initiate the procedure.

• When gelatinous pellet is detected in the semen, leave it at room temperature for 30 more minutes or remove those using filters or other appropriate equipment.

• Warm SepaSperm® Solution to 37°C prior to use. The products should be tightly capped when used in CO₂ incubator to prevent pH alteration.

• Properly mix the prepared density gradients by inverting the bottles several times.

Instructions for use:

1. Prepare a dual gradient system (45% – 90% or 40% – 80%) for each semen sample with a rinsing medium.

• 90% gradient: Dilute SepaSperm® Solution with 10% of a rinsing medium.

• 45% gradient: Dilute SepaSperm® Solution with 55% of a rinsing medium.

NOTE: Gradients should be prepared and repacked under sterile conditions (e.g. LAF bench ISO class 5 and sterile recipient). For optimal results:

• Prepare the gradients maximum 24 hours prior to use and store at 2-8°C and warm the gradients to room temperature or 37°C one hour before use.

• Prepare the gradients needed for one day in one time and mix well. Discard excess unused prepared gradient.

2. Carefully place 2.5 ml of 90% density gradient on the bottom of the centrifuge tube.

3. Add 2.5 ml of 45% density gradient on top of the 90% density gradient.

4. Gently place the liquefied semen on the top layer (45% density gradient).

5. Centrifuge the tubes at 350 to 400 g for 15 to 18 minutes. If the pellets are not visible, centrifuge for a further 3 minutes.

6. Remove the supernatant.

QUALITY CONTROL SPECIFICATIONS:

The following tests were performed for each lot of this product:

- Sterility: SAL 10⁻³
- Endotoxin: < 0.5 EU/ml
- Sperm Survival test: ≥ 80% survival after 4 hours exposure; ≥ 75% survival after 24 hours exposure.
- pH: 7.20 – 7.90 (Lot release criteria: 7.20-7.60)
- Osmolality: 300-330 mOsm/kg

All results are reported on Certificate of Analysis which is available upon request.

STORAGE INSTRUCTIONS AND**STABILITY:**

- Store between 2-8°C.
- Keep away from sunlight.
- Once opened, store the product in the original packaging.
- Once opened, the product is stable for at least 7 days under aseptic conditions and store between 2-8°C. Discard the remaining media if this period has passed.

• The product is stable after transport (max. 5 days) at elevated temperatures (< 37°C).

• This product is stable until the expiration date labelled on the bottle if the packaging is unopened or undamaged

CONTRAINDICATIONS:

- Do not re-sterilize.
- Do not freeze the product.
- Do not use the product after expiration date.
- Do not use if packaging is damaged or broken.
- Do not use if product becomes discolored, cloudy or shows any evidence of microbial contamination.

WARNINGS:

• Always works under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid contamination. Aseptic technique should be used

• In case of eye or skin contact with SepaSperm® Solution, immediately flush eye/skin with water.

• Observe all federal, state and local environmental regulations when discarding the product.

• In case of infection, dispose the product appropriately in a prescribed manner.

• The user shall be responsible for any problems caused by non-conformity to the present IFU.

• This product is intended to be used by medical specialist trained in fertility treatment.

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