

Results of Cryotop SC - Clinical Trial

	Cryotop SC System	Cryotop Open System
Nº Randomized donors	185	185
Nº Donations (Recipients)	159	146
Oocytes received (mean ± SD)	2,048	1,787
Survival rate	1,787 (87.3%)	1,639 (91.7%)
Fertilization rate (2PN)	1387 (74.6%)	1220 (74.4%)
Top quality blastocyst	196/561 (34.9%)	201/571 (35.2%)
Cycles with/cryopreserv. surplus embryos	125 (78.6%)	118 (74.7%)
Nº transfers	151 (94.9%)	134 (91.7%)
Clinical pregnancy/cycle	89 (56.0%)	68 (46.6%)